

# Blotting Protocol

- 1] **Soak** gel in 0.5M NaOH for 30min on rocker
- 2] **Soak** gel in 10X SSC for 15min on rocker
- 3] **Cut** one sheet of nylon membrane and two sheets of Whatmann 3MM paper to the size of the gel
- 4] **Cut** two lengths of whatmann paper wider than the gel and long enough to fit from end to end in transfer buffer tray
- 5] **Wet** the nylon membrane in water
- 6] **Soak** the nylon membrane and whatmann paper in 20X SSC for 1-2min
- 7] **Set** up the capillary transfer
  - fill large pyrex dish (cleaned with RNaseZap) with 20X SSC
  - place a support in the middle
  - place the two long whatmann wicks over the support and remove air bubbles
  - place gel upside down on whatmann wick
  - place plastic wrap around the edges of gel to prevent short circuits
  - place presoaked nylon membrane over gel and remove air bubbles
  - place the two presoaked sheets of whatmann over nylon membrane and remove air bubbles
  - place a stack of 15-20cm of dry/cut to gel size paper towels
  - place a piece of glass on top of paper towels
  - place 1kg weight on top of glass
  - let proceed for 12-18hr, replacing the wet paper towels at least once
- 8] **Disassemble** the capillary transfer set up
- 9] **Carefully** peel away the gel from the membrane
- 10] **Mark** the lanes with a blunt ended pencil
- 11] **Check** transfer on uv box and mark in pencil the ladder positions
- 12] **Soak** membrane in 3X SSC for 10min
- 13] **Autocrosslink** with stratalinker in Proctor Lab
- 14] **Place** in seal a meal bag if not using right away and store at 4C

