

## Radioactive Probe Preparation

**Order** Prime-a-Gene Labeling System from Promega (Cat # U1100)

**Order** from CORD Perkin Elmer Life Science alpha-32P-dCTP, 3000 CI/mmol, 10mCI/ml (Cat # BLU-513H)

- 1] **PCR** the gene of interest and clean up with Qiagen PCR purification kit
- 2] **Spec** to determine concentration (usually 10-50ng/ul)
- 3] **Place 25ng** of clean PCR rxn in 0.5ml tube and bring up to 5ul with nuclease free H<sub>2</sub>O
- 4] **Denature** for 2 min at 95C in PCR Machine
- 5] **Rapidly chill** in ice water bath for at least 2 min
- 6] **Spin** down to bring everything to bottom
- 7] **Make dNTP mixture**
  - from the kit, mix an equal amount of dATP, dGTP, dTTP
- 8] **Assemble reaction in the order shown in a 0.5ml tube**
  - 25ul Nuclease free H<sub>2</sub>O
  - 10ul 5X Labeling buffer
  - 2ul of dNTP mixture (dATP, dGTP, dTTP)
  - 5ul of denatured template
  - 2ul Nuclease free BSA
  - 5ul alpha-32P-dCTP
  - 1ul Klenow Polymerase
- 9] **Mix** gently and centrifuge briefly to bring down contents
- 10] **Incubate** reaction at room temp for 60min
- 11] **Terminate** reaction by heating to 95C for 2min and rapidly chilling in ice water bath
- 12] **Add** 2.083ul of 0.5M EDTA
- 13] **Prepare G-50** spin column by spinning in 1.5ml tube at 2Xg for 1 min
- 14] **Blot** the extra drop of water off of the end of the column after spinning and place in a new 1.5ml tube
- 15] **Load** the 50ul probe sample above directly to the middle of the slanted G50 column
- 16] **Spin** at 2Xg for 1 min
- 17] **Take one microliter** of sample and check with Geiger counter