

XTT Assay

Day 1:

Start 5 ml YPD culture
Set up template in notebook
Check to be sure reagents are available

Day 2:

Perform hemocytometer counts: cells/5 small boxes x 5×10^4
Dilute cells to 10^6 in RPMI (record bottle used in notes)
Place 100 μ l in each well
-Include controls without cells
-Do not use outside wells, except for controls
Place 200 μ l saline in outside wells
Parafilm and place in foil
Incubate 37° for 24 hours

Day 3:

Prepare antifungal to have ready for after washes
Remove media
-If using same tips, must have same cell type in each well.
-Do not contaminate controls.
Wash gently with 100 μ l PBS 2-3 x
-Be sure to remove all fluid each time
-New tip for each wash
-Place tip in corner and be careful not too disrupt the biofilm
Add 90 or 100 μ l antifungal
-If making serial dilutions, must dilute in another plate to avoid transfer of cells.
-If making serial dilutions in another plate, transfer only 90 μ l to biofilms.
Add 90 or 100 μ l RPMI
-Add only 90 μ l of RPMI if 90 μ l of antifungal added
Incubate 24 hours at 37°

Day 4

Thaw XTT/PMS or FDA
-Place in warm water if needed quickly
Centrifuge to remove sediment
Remove media
-Use new tip each well.
Wash gently with 100 μ l PBS 2-3 x
-Be sure to remove all fluid each time
-New tip for each wash
-Place tip in corner and be careful not too disrupt the biofilm
Turn off direct lights
Add 90 μ l XTT (0.5mg/ml or 1 mg/ml)
-Pipette each time, do not use multichannel
Add 10 μ l PMS
-Okay to use multichannel, but only need to add about 6ml to Petri dish.
Place in foil and Incubate at 37° 2 hours
Read at 490 nm

XTT (1mg/ml):

1mg/ml 1x PBS

Centrifuge to remove sediment

Collect supernatant

Store -20° in foil

PMS (0.32 mg/ ml):

Add 9.6mg PMS to 30 ml sterile water

Place directly into foil and store -20°