

**XTT drug susceptibility assay for *C. albicans* biofilms formed in 96-well plates (flat-bottom wells; Costar 3590/3361) [updated 4-3-2019 by Hiram Sanchez]**

-Culture on warmed YPD (yeast-peptone-dextrose) plates (3-phase streak) from -80°C glycerol stocks- incubate at **30°C for 1-2 days** (after this keep plates in 4°C for 1-2 weeks)

Day 1:

1. From a single colony grow your samples overnight in 15mL conical tubes containing 5 ml of YPD, **incubate at 30°C with gently shaking 50 to 200 rpm**. Include reference strain (i.e. Wild Type-SN250).
  - 4 samples fit per plate, 1 reference strain and 3 tested strains

Day 2:

1. Remove cultures and pulse vortex for proper mixing.
2. To count cells with hemocytometer see addendum for details.
3. Using a 96 well flat-bottom plate (Costar 3590/3361) label and inoculate with the plate with **1 X10<sup>6</sup> cell/ml** as described (see below figure)
  - Place **200µL of 1x PBS** in Rows A and H and column 12 to prevent evaporation of samples
  - Place **100µL of each strain** to be tested into their corresponding quadrant
  - Place **100µL of RPMI** in column 1 as a negative control

Plate Setup (Fluconazole)

	1	2	3	4	5	6	7	8	9	10	11	12
A	PBS 200µL											
B	RPMI 100µL	Wild type 100µL					Sample B 100µL					PBS 200µL
C												
D												
E		Sample A 100µL					Sample C 100µL					
F												
G												
H	PBS 200µL											

4. Seal plate with parafilm. Label with your initials, time/date, and testing drug.
5. Place plate in **37°C incubator**, without shaking for **6h**.
6. After **6hrs** of incubation remove plate from incubator, slowly remove media and dose with fluconazole-try not to disrupt biofilm. **(DOSE 1)**
7. Seal plate with parafilm and **incubate at 37°C for 24hrs**.

### Day 3:

1. **24hrs** after the first fluconazole dose-slowly remove media (**put in 96-well plate for dispersion**) and dose all wells with fluconazole a second time-try not to disrupt biofilm. (**DOSE 2**)
2. Parafilm; **incubate at 37°C without shaking for 24 hr.**

### Day 4:

1. Slowly remove media (**put in 96-well plate for dispersion**) and wash each well carefully (unless fragile biofilm-then don't wash), so as not to disturb the biofilm, with **100µL of 1x PBS one time.**
2. Combine **1mL PMS+9mL XTT solution** (aliquots kept in -20°C-both are light sensitive)-**only enough for one plate.**  
Add **100µL of XTT+PMS** solution to each well.
3. Wrap plate in foil and incubate at **37°C without shaking for 30 min.**
4. Read absorbance at  $\lambda$ : 490 nm-**4 sample Fluconazole XTT.prt** (protocol to use with Gen5)

### Dosing

1. Slowly remove media from all wells.
2. Add **180µL of drug+RPMI** (1:2 dilution-that's why final concentration is half of stock) to wells slowly-try not to disrupt biofilm. (**if using same tips plate lowest concentration to highest**)
3. When Finished parafilm plate and place in **37°C incubator.**

### **Reagents:**

**XTT:** (Sigma X4626; kept at 4°C).

Prepare at 0.83 mg/mL in 1x PBS.

Add 7.5 mg of XTT salt to 9mL 1x PBS

Centrifuge (3000 RPM for 1 min at RT) or filter to remove sediment

Collect supernatant

Foil tube and store at -20°C

**Phenazine methosulfate (PMS):** (Sigma P9625; kept at -20°C)

Prepare at 0.32 mg/mL in dH<sub>2</sub>O.

Add 3.2 mg of PMS to 10mL ddH<sub>2</sub>O (filter-sterilize)

Shake to dissolve and place in foil and store at -20°C.

**RPMI + MOPS:** (RPMI-1640 Medium HyClone SH30011.04; kept at 4°C, MOPS BP308; kept at RT)

Add 10.4 g RPMI to 700 mL of ddH<sub>2</sub>O into a 1-2L beaker with a magnetic stir bar. Add 34.5 g of MOPS and mix well (with magnetic stir bar). Adjust pH to 7.0 using 5N liquid NaOH or NaOH pellets. Bring volume to 1 L and filter sterilize using Polyethersulfone Membrane (PES) 0.22 µm pore size filters into 500 mL bottles. Store at 37°C or RT.

### Plate Layout (Fluconazole)


	1	2	3	4	5	6	7	8	9	10	11	12		
A	PBS 200µL													
B	RPMI 100µL	Ctrl Wells (No Drug)	Highest Drug Concentration	Wild Type 100µL			Ctrl Wells (No Drug)	Highest Drug Concentration	Sample B 100µL			PBS 200µL		
C														
D														
E					Sample A 100µL								Sample C 100µL	
F					Decreasing Drug Concentration								Decreasing Drug Concentration	
G														
H			PBS 200µL											

### Dosing Plate (Fluconazole)










	1	2	3	4	5	6	7	8	9	10	11	12	
A			1000µg/ml F.C.	500	250	125		1000	500	250	125		
B	RPMI 100µL	Fluconazole 200µL (2 mg/mL) (2000µg/mL)		RPMI 100µL			Fluconazole 200µL (2 mg/mL)		RPMI 100µL				
C													
D													
E													
F													
G													
H						100µL			100µL	100µL			100µL

**+100 µL RPMI before plating**

**Plate Layout (Amphotericin, Micafungin, Biocides, Etc.)**

	1	2	3	4	5	6	7	8	9	10	11	12	
A	PBS 200µL												
B	RPMI 100µL	Ctrl Wells (No Drug)	Highest Drug Concentration	Wild Type 100µL						PBS 200µL			
C													
D													
E				Sample A 100µL									
F													
G				Decreasing Drug Concentration 									
H	PBS 200µL												

**Dosing Plate (Amphotericin, Micafungin, Biocides, Etc.)**

	1	2	3	4	5	6	7	8	9	10	11	12	
A													
B	Saline 100µL		Drug Being Tested 200µL (Highest Concentration)	Saline 100µL									
C													
D													
E													
F													
G													
H													

Discard 100µL

# Addendum

## Cell counting with hemocytometer

- 1- Prepare 5 ml of a 1:100 dilution of your overnight culture (50  $\mu$ l of overnight culture into 4,950  $\mu$ l of RPMI Media) mix well.
- 2- Take 10  $\mu$ l of the cell dilution and place it in the hemocytometer.
- 3- Count the number of cells in 4 large squares (you can also count the cells in 4 medium and 4 small size squares but all you need is 4 of the same size, you decide what to use) see Fig 1.
- 4- Average the number of cells in all of the four squares and multiply the average by 1000  $\mu$ l.
- 5- Then divide that number by one of these three;
  - a. If large squares were used divide by 0.00625 mm<sup>3</sup>
  - b. If medium squares divide by 0.004 mm<sup>3</sup>
  - c. If small squares divide by 0.00025 mm<sup>3</sup>

For example:

Dilution was 1:100 and you counted the cells in four large squares, in each individual square you counted 15,19,18 & 17, when added and divided by 4, you end up with  $69/4 = 17.25$   
17.25 is multiply by 1000  $\mu$ l = 17,250. Since you counted the cells in the large square you need to divide the 17,250 by 0.00625.  $17,250/0.00625 = 2760000$  which is the same as  $2.76 \times 10^6$  cell/ ml in your 1:100 dilution. The total number of cells in your original overnight culture is  $2.76 \times 10^6/ \text{ml} \times 100 = 2.76 \times 10^8$  cells/ml.

So for your experiment you need **1 X10<sup>6</sup> cell/ml**, so follow this;

(Number of cells you want) x (the volume in ml you need)

Total number of cells you have in your (1:100 dilution)

$$= (1 \times 10^6) \times (5 \text{ ml}) / 2.76 \times 10^6 \text{ per ml} = 1.81 \text{ ml}$$

Add 1.81 ml of your 1:100 dilution to 3.18 ml of media to make 5 ml at 1 X10<sup>6</sup> cell/ml.

Fig.1

